

## AUTOCLAV.DOUBLE AD CONICAL ARM CUFF(47X16CM) W/CONNECTOR F/EGP200

Réf. GHP552813

Hpm™



# AUTOCLAV.DOUBLE AD CONICAL ARM CUFF(47X16CM) W/CONNECTOR F/EGP200



Réf. **GHP552813**

Nylon cuff with double silicone bladder suitable for autoclaving. With male Luer lock connector and cap for cleaning and sterilization.

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3700446033774
Guarantee (month):	12	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,615
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:05 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference GHP552813