

## ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC

Réf. 527500



CE  
0459



# ARM BLOOD PRESSURE MONITOR SPENGLER TENSONIC



Réf. 527500

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe IIa              |               |
| Marquage                   | CE0459    | Barcode                 | 3700446022624 |
| Guarantee (month):         | 36        | Customs Code            | 90189010      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              | 0,765     | Weight packaging (kg) : | 0,765         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | -20       | Temperature max         | 55            |
| Humidity min               | 15        | Humidity max            | 90            |



08/12/2022 12:37 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 527500

All characteristics are given for information only and may be changed