

## DISPOSABLE SPECULUM DIA. 4 MM, BOX OF 250, BLACK, ADULT, SPENGLER

Réf. L20748



# DISPOSABLE SPECULUM DIA. 4 MM, BOX OF 250, BLACK, ADULT, SPENGLER



Réf. **L20748**

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | Yes       | Classe I                |               |
| Marquage                   | CE        | Barcode                 | 3700446031275 |
| Guarantee (month):         | 6         | Customs Code            | 90330000      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   | 9,6           |
| Length (cm) :              |           | Length packaging (cm) : | 11,7          |
| Height (cm) :              |           | Height packaging (cm) : | 11,4          |
| Weight (kg) :              |           | Weight packaging (kg) : | 0,135         |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



31/03/2022 01:24 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference L20748

All characteristics are given for information only and may be changed