

GLASS SPENGLER (DIA.167MM) MAXI+2



Réf. 310270



GLASS SPENGLER (DIA.167MM) MAXI+2



Réf. 310270

| FEATURES | | | |
|----------------------------|-----------|-------------------------|---------------|
| Medical device | No | | |
| Marquage | CE | Barcode | 3700446025076 |
| Guarantee (month): | 60 | Customs Code | 9018901000 |
| Sterile | No | | |
| INFORMATIONS | Main unit | Packaged unit | |
| Width (cm) : | | Width packaging (cm): | 18,0 |
| Length (cm) : | | Length packaging (cm) : | 18,0 |
| Height (cm) : | | Height packaging (cm) : | 2,0 |
| Weight (kg) : | | Weight packaging (kg) : | 0,086 |
| SPECIFIC STORAGE CONDITION | | | |
| Temperature min | (N/A) | Temperature max | (N/A) |
| Humidity min | (N/A) | Humidity max | (N/A) |



31/03/2022 02:21 - This article is available on our website : www.gsh-med.fr with the following reference 310270

All characteristics are given for information only and may be changed