

## CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/LARGE ADULT (L)

Réf. 337850



# CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/LARGE ADULT (L)



Réf. **337850**

Single tube Clinic® cuff for Lian® + clamp compatible w/ Nihon Kohden patient monitors

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446030872
Guarantee (month):	6		Customs Code	9018901000
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :		Width packaging (cm) :		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :		
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



26/04/2022 12:33 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 337850

All characteristics are given for information only and may be changed