

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/CHILD (S)



Réf. 337855



CLINIC CUFF/1 TUBE/WITH COUPLING # 316 152 AND CLAMP/CHILD (S)



Réf. **337855**

Single tube Clinic® cuff for Lian® + clamp compatible w/ Nihon Kohden patient monitors

FEATURES			
Medical device	No		
Marquage		Barcode	3700446030889
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference 337855

All characteristics are given for information only and may be changed