

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 163 AND CLAMP/LARGE ADULT (XL)

Réf. 337881



CLINIC CUFF/1 TUBE/WITH COUPLING # 316 163 AND CLAMP/LARGE ADULT (XL)



Réf. 337881

Single tube Clinic® cuff for Lian® + clamp compatible w/ Datascope monitor

FEATURES			
Medical device	Yes	Classe I	
Marquage			
Guarantee (month):	6	Barcode	3700446030964
		Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference 337881

All characteristics are given for information only and may be changed