

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 163 AND CLAMP/ADULT (M)

Réf. 337883



CLINIC CUFF/1 TUBE/WITH COUPLING # 316 163 AND CLAMP/ADULT (M)



Réf. 337883

Single tube Clinic® cuff for Lian® + clamp compatible w/ Datascope monitor

FEATURES			
Medical device	Yes	Classe I	
Marquage		Barcode	3700446000172
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit		Packaged unit
Width (cm) :			Width packaging (cm):
Length (cm) :			Length packaging (cm) :
Height (cm) :			Height packaging (cm) :
Weight (kg) :			Weight packaging (kg) : 0,125
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference 337883

All characteristics are given for information only and may be changed