

NIPB TUBING 2,5M FOR 4VS MONITOR



Réf. 551624



NIPB TUBING 2,5M FOR 4VS MONITOR



Réf. 551624

FEATURES			
Medical device	No		
Marquage		Barcode	3661040213842
Guarantee (month):	6	Customs Code	85049099
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	16,0
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,118
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:31 - This article is available on our website : www.gsh-med.fr with the following reference 551624

All characteristics are given for information only and may be changed