

DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSONIC

Réf. 527500P



DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSORIC



Réf. 527500P

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3661040213774
Guarantee (month):	36	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min	15	Humidity max	90



31/03/2022 02:37 - This article is available on our website : www.gsh-med.fr with the following reference 527500P

DISPLAY OF 6 ARM BLOOD PRESSURE MONITOR SPENGLER TENSONIC



Réf. 527500P

COMPOSITION

Code	Descriptive	Quantity
527500	Arm blood pressure monitor Spengler Tensonic	6