

BLACK NYLON SOFT BAG SPENGLER, MEDIUM

Réf. 535020



BLACK NYLON SOFT BAG SPENGLER, MEDIUM



Réf. 535020

Nylon pouch 18x13x6 cm with éclair® zip closure

FEATURES			
Medical device	No		
Marquage		Barcode	3700446032395
Guarantee (month):	60	Customs Code	90330000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,085
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:50 - This article is available on our website : www.gsh-med.fr with the following reference 535020

All characteristics are given for information only and may be changed