

SPO2 REUSABLE DIGITAL SENSOR FOR ADULT

Réf. **DSU250900**

| FEATURES | | | | | |
|----------------------------|--------|-------------------------|----------------------|-----------------|-------|
| Medical device | Yes | Classe I | Accuracy measurement | Classe I | |
| Marquage | CE0413 | | Barcode | 010869887032265 | |
| Guarantee (month): | 0 | | Customs Code | 90181990 | |
| Sterile | No | | | | |
| INFORMATIONS | | Main unit | Packaged unit | | |
| Width (cm) : | | Width packaging (cm): | | | |
| Length (cm) : | | Length packaging (cm) : | | | |
| Height (cm) : | | Height packaging (cm) : | | | |
| Weight (kg) : | | Weight packaging (kg) : | | 0,052 | |
| SPECIFIC STORAGE CONDITION | | | | | |
| Temperature min | | (N/A) | Temperature max | | (N/A) |
| Humidity min | | (N/A) | Humidity max | | (N/A) |



15/11/2022 12:46 - This article is available on our website : www.gsh-med.fr with the following reference DSU250900

All characteristics are given for information only and may be changed