

Réf. DSU551463



Réf. DSU551463

FEATURES				
Medical device	Yes	Classe IIa	Accuracy measurement	Classe II
Marquage	CE0413		Barcode	0120884521008134
Guarantee (month):	3		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,0	
Length (cm) :		Length packaging (cm) :	15,0	
Height (cm) :		Height packaging (cm) :	5,0	
Weight (kg) :		Weight packaging (kg) :	0,086	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 02:03 - This article is available on our website : www.gsh-med.fr with the following reference DSU551463

All characteristics are given for information only and may be changed