

## BLACK NYLON SOFT BAG SPENGLER, LARGE

Réf. 535019



# BLACK NYLON SOFT BAG SPENGLER, LARGE



**Réf. 535019**

Nylon pouch 25x20x5 cm with éclair® zip closure

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446002008
Guarantee (month):	60		Customs Code	90330000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,111	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 01:21 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 535019

All characteristics are given for information only and may be changed