

SMARTLIGHT OTOSCOPE BLACK

Réf. 570315



SMARTLIGHT OTOSCOPE BLACK



Réf. 570315

Conventional lighting. Lightweight and robust, this otoscope is perfect for ENT auscultations especially on the move. The XHL lamp diffuses a concentrated and homogenous light beam (optional LED lamp available). Fit with an ABS and Brass head, SMARTLIGHT also comes with end clamp and built-in switch, LR6 AA batteries included.

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE	Barcode	3700446041557
Guarantee (month):	12	Customs Code	90189042
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	19,0
Length (cm) :		Length packaging (cm) :	10,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,230
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	10	Humidity max	90



31/03/2022 02:46 - This article is available on our website : www.gsh-med.fr with the following reference 570315

All characteristics are given for information only and may be changed