

TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER

Réf. SP3105



TEMP DUO II SPENGLER FOREHEAD AND AURICULAR THERMOMETER



Réf. **SP3105**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446041786
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,5
Length (cm) :		Length packaging (cm) :	18,5
Height (cm) :		Height packaging (cm) :	5,0
Weight (kg) :		Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:52 - This article is available on our website : www.gsh-med.fr with the following reference SP3105

All characteristics are given for information only and may be changed