

HANDLE FOR SMARTLED 5500-R



Réf. 579531

HANDLE FOR SMARTLED 5500-R



Réf. 579531

FEATURES			
Medical device	No		
Marquage		Barcode	3661040213910
Guarantee (month):	0	Customs Code	90189090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	2,0
Length (cm) :		Length packaging (cm) :	12,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	0,050
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:35 - This article is available on our website : www.gsh-med.fr with the following reference 579531

All characteristics are given for information only and may be changed