

LED BULB FOR LITESTICK



Réf. 570690

LED BULB FOR LITESTICK



Réf. 570690

Spare part for Litestick® penlight.

FEATURES			
Medical device	No		
Marquage		Barcode	3700446009908
Guarantee (month):	0	Customs Code	90189084
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/08/2022 12:43 - This article is available on our website : www.gsh-med.fr with the following reference 570690

All characteristics are given for information only and may be changed