

GLASS FOR SPENGLER VAQUEZ LAUBRY COLOUR/HOPITAL (LARGE DIAL) (X12)



Réf. 310220



GLASS FOR SPENGLER VAQUEZ LAUBRY COLOUR/HOPITAL (LARGE DIAL) (X12)



Réf. 310220

FEATURES			
Medical device	No		
Marquage		Barcode	3700446024901
Guarantee (month):	60	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	2,0
Length (cm) :		Length packaging (cm) :	2,0
Height (cm) :		Height packaging (cm) :	2,0
Weight (kg) :		Weight packaging (kg) :	1,260
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:21 - This article is available on our website : www.gsh-med.fr with the following reference 310220

All characteristics are given for information only and may be changed