

## LIAN CUFF/1 TUBE/NYLON/BLACK OXFORD/CHILD (S)

Réf. 333020



# LIAN CUFF/1 TUBE/NYLON/BLACK OXFORD/CHILD (S)



**Réf. 333020**

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color.  
Washable at 30-40° C.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446008741
Guarantee (month):	6		Customs Code	9018901000
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm) :	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,100
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 02:23 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 333020

All characteristics are given for information only and may be changed