

Réf. DSU230723



Réf. DSU230723

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413		
Guarantee (month):	24	Accuracy measurement	Classe I
		Barcode	3700446012328
Sterile	No	Customs Code	
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,101
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU230723

All characteristics are given for information only and may be changed