

RECHARGING DOCK FOR SMARTLED 5500-R



Réf. 579534

RECHARGING DOCK FOR SMARTLED 5500-R



Réf. 579534

FEATURES			
Medical device	No		
Marquage		Barcode	3700446032913
Guarantee (month):	0	Customs Code	
Sterile	No		
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	14,0
Length (cm) :		Length packaging (cm) :	17,0
Height (cm) :		Height packaging (cm) :	8,0
Weight (kg) :	0,354	Weight packaging (kg) :	0,354
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:55 - This article is available on our website : www.gsh-med.fr with the following reference 579534

All characteristics are given for information only and may be changed