

## REPLACEMENT HEAD BLUEBERRY FOR OTOSCOPE SMARTLED 5500



Réf. 570518B

# REPLACEMENT HEAD BLUEBERRY FOR OTOSCOPE SMARTLED 5500



Réf. **570518B**

| FEATURES                   |           |                         |               |
|----------------------------|-----------|-------------------------|---------------|
| Medical device             | No        |                         |               |
| Marquage                   |           | Barcode                 | 3700446041779 |
| Guarantee (month):         | 0         | Customs Code            | 90189042      |
| Sterile                    | No        |                         |               |
| INFORMATIONS               | Main unit | Packaged unit           |               |
| Width (cm) :               |           | Width packaging (cm):   |               |
| Length (cm) :              |           | Length packaging (cm) : |               |
| Height (cm) :              |           | Height packaging (cm) : |               |
| Weight (kg) :              |           | Weight packaging (kg) : |               |
| SPECIFIC STORAGE CONDITION |           |                         |               |
| Temperature min            | (N/A)     | Temperature max         | (N/A)         |
| Humidity min               | (N/A)     | Humidity max            | (N/A)         |



31/03/2022 02:56 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 570518B

All characteristics are given for information only and may be changed