

REPLACEMENT HEAD SOFT GREEN FOR OTOSCOPE SMARTLED 5500 PEDIA



Réf. 570518P

REPLACEMENT HEAD SOFT GREEN FOR OTOSCOPE SMARTLED 5500 PEDIA



Réf. 570518P

FEATURES			
Medical device	No		
Marquage		Barcode	3700446041755
Guarantee (month):	0	Customs Code	90189042
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:56 - This article is available on our website : www.gsh-med.fr with the following reference 570518P

All characteristics are given for information only and may be changed