

BABY TO TODDLER ELECTRONIC SCALE



Réf. EHANDI3



BABY TO TODDLER ELECTRONIC SCALE



Réf. **EHANDI3**

FEATURES			
Medical device	No		
Marquage	CE	Barcode	3661040216669
Guarantee (month):	12	Customs Code	84231090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	34,8
Length (cm) :		Length packaging (cm) :	63,5
Height (cm) :		Height packaging (cm) :	16,9
Weight (kg) :		Weight packaging (kg) :	1,500
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:53 - This article is available on our website : www.gsh-med.fr with the following reference EHANDI3

All characteristics are given for information only and may be changed