

LIAN METAL HAND-HELD MANOMETER

Réf. 310045



LIAN METAL HAND-HELD MANOMETER



Réf. 310045

Dials 0 to 300 mmHg

FEATURES				
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement	Classe I
Marquage			Barcode	3700446000141
Guarantee (month): 60			Customs Code	90189010
Sterile	No			
INFORMATIONS	Main unit		Packaged unit	
Width (cm) :			Width packaging (cm):	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,380
SPECIFIC STORAGE CONDITION				
Temperature min		-20	Temperature max	70
Humidity min		20	Humidity max	85



31/03/2022 02:21 - This article is available on our website : www.gsh-med.fr with the following reference 310045

All characteristics are given for information only and may be changed