

LIAN NANO HAND-HELD MANOMETER/BLACK CASE/SILVER RING

Réf. 310065



LIAN NANO HAND-HELD MANOMETER/BLACK CASE/SILVER RING



Réf. 310065

Dials 0 to 300 mmHg

FEATURES			
Medical device	Yes	Classe Im (fonction de mesurage)	Accuracy measurement Classe I
Marquage	CE0459	Barcode	3700446023805
Guarantee (month):	60	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,423
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	70
Humidity min	20	Humidity max	85



31/03/2022 02:21 - This article is available on our website : www.gsh-med.fr with the following reference 310065

All characteristics are given for information only and may be changed