

OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY

Réf. 221052



OXYSTART® - DIGITAL FINGERTIP PULSE OXIMETER BLUBERRY



Réf. **221052**

FEATURES			
Medical device	Yes	Classe IIb	
Marquage	CE0123	Barcode	3700446041809
Guarantee (month):	24	Customs Code	90181910
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,135
SPECIFIC STORAGE CONDITION			
Temperature min	-40	Temperature max	60
Humidity min	5	Humidity max	95



31/05/2022 12:31 - This article is available on our website : www.gsh-med.fr with the following reference 221052

All characteristics are given for information only and may be changed