

6 X BLACK BULBS WITH VALVE FOR VL CLASSIC, VL CLINIC, VL NANO, MAXI+3

Réf. 311014



6 X BLACK BULBS WITH VALVE FOR VL CLASSIC, VL CLINIC, VL NANO, MAXI+3



Réf. **311014**

Classic® Black bulbs/Classic® valves

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446026950
Guarantee (month):	60		Customs Code	9018901000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,500	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



23/06/2022 21:04 - This article is available on our website : www.gsh-med.fr with the following reference 311014

All characteristics are given for information only and may be changed