

AC/DC ADAPTOR 7.5 VOLTS FOR ES 60



Réf. 550615



CE
0044

AC/DC ADAPTOR 7.5 VOLTS FOR ES 60



Réf. 550615

FEATURES			
Medical device	No		
Marquage	CE0044	Barcode	3700446008918
Guarantee (month):	12	Customs Code	8504318090
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 01:52 - This article is available on our website : www.gsh-med.fr with the following reference 550615

All characteristics are given for information only and may be changed