

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/CHILD (S)

Réf. 333029



LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/CHILD (S)



Réf. 333029

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446006877
Guarantee (month):	6		Customs Code	9018901000
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm) :		
Length (cm) :		Length packaging (cm) :		
Height (cm) :		Height packaging (cm) :		
Weight (kg) :		Weight packaging (kg) :	0,081	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



25/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 333029

All characteristics are given for information only and may be changed