

LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/INFANT (XS)

Réf. 333039



LIAN CUFF/1 TUBE/COTTON/CARTOON DESIGN/INFANT (XS)



Réf. 333039

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

| FEATURES | | | | |
|----------------------------|-----------|-------------------------|----------------------|---------------|
| Medical device | Yes | Classe I | Accuracy measurement | Classe I |
| Marquage | | | Barcode | 3700446000103 |
| Guarantee (month): | 6 | | Customs Code | 9018901000 |
| Sterile | No | | | |
| INFORMATIONS | Main unit | Packaged unit | | |
| Width (cm) : | | Width packaging (cm) : | | |
| Length (cm) : | | Length packaging (cm) : | | |
| Height (cm) : | | Height packaging (cm) : | | |
| Weight (kg) : | | Weight packaging (kg) : | 0,100 | |
| SPECIFIC STORAGE CONDITION | | | | |
| Temperature min | (N/A) | Temperature max | (N/A) | |
| Humidity min | (N/A) | Humidity max | (N/A) | |



25/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference 333039

All characteristics are given for information only and may be changed