

CLINIC CUFF/1 TUBE/WITH COUPLING # 316 162 AND CLAMP/ADULT (M)

Réf. 337862



CLINIC CUFF/1 TUBE/WITH COUPLING # 316 162 AND CLAMP/ADULT (M)



Réf. **337862**

Single tube Clinic® cuff for Lian® + clamp compatible w/ Datascope, Drager, HP/ Agilent/ Philips, Marquette, Hellige, Nellcor, Siemens, Spacelabs patient

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446001940
Guarantee (month):	6		Customs Code	9018901000
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm) :	
Length (cm) :			Length packaging (cm) :	
Height (cm) :			Height packaging (cm) :	
Weight (kg) :			Weight packaging (kg) :	0,130
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



26/04/2022 12:34 - This article is available on our website : www.gsh-med.fr with the following reference 337862

All characteristics are given for information only and may be changed