

FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

Réf. DSU250904

 SunTech Medical®

CE
0413



FLEXIBLE CHILD DIGITAL SPO2 SENSOR (90CM CABLE)

Réf. **DSU250904**

Flexible child digital SpO2 sensor (90 cm cable)

FEATURES			
Medical device	Yes	Classe I	Accuracy measurement
Marquage	CE0413		Classe I
Guarantee (month):	0	Barcode	8698870322296
Sterile	No	Customs Code	90181990
INFORMATIONS		Main unit	
		Packaged unit	
Width (cm) :		Width packaging (cm):	2,0
Length (cm) :		Length packaging (cm) :	14,0
Height (cm) :		Height packaging (cm) :	16,0
Weight (kg) :		Weight packaging (kg) :	0,010
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:02 - This article is available on our website : www.gsh-med.fr with the following reference DSU250904

All characteristics are given for information only and may be changed