

## GREY SPENGLER (XL) CUFF (32-48CM) FOR TENSONIC & AUTOTENSIO

Réf. 527465



# GREY SPENGLER (XL) CUFF (32-48CM) FOR TENSORIC & AUTOTENSIO



Réf. **527465**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0197	Barcode	3700446041915
Guarantee (month):	6	Customs Code	90262040
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	13,0
Length (cm) :		Length packaging (cm) :	16,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,130
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



26/04/2022 12:35 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference 527465

All characteristics are given for information only and may be changed