

GREY SPENGLER (ML) CUFF (22-42CM) FOR TENSORIC & AUTOTENSIO

Réf. 527464



GREY SPENGLER (ML) CUFF (22-42CM) FOR TENSONIC & AUTOTENSIO



Réf. 527464

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0197	Barcode	3700446041908
Guarantee (month):	24	Customs Code	90262040
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,102
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:01 - This article is available on our website : www.gsh-med.fr with the following reference 527464

All characteristics are given for information only and may be changed