

NON-CONTACT THERMOMETER TEMPO PRO SPENGLER

Réf. SP3110



NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Réf. SP3110

Infrared (IR), non-contact, clinical thermometer, digital, handheld, range 34 # 43 °C

At a distance of a few centimetres from the temple, its display instantly shows the temperature.
Without contact to avoid cross-contamination, Tempo Pro is easy to disinfect and always stays clean.

TECHNICAL CHARACTERISTICS

- Measurement : Infrared
- Distance de mesure : - 5 cm
- Unit of measurement : Celsius or Fahrenheit
- Measuring range : 32,0°C – 43,0°C / 89,6°F - 109,4°F
- Resolution : 0,1°C / F°
- Forehead accuracy : 0,2°C / 0,4F°
- Memories : 30 memories
- Measurement time: Less than 1 second
- Screen : Backlit LCD screen
- Auto power off
- 2 x 1,5V AAA batteries (included)

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0123	Barcode	3700446041823
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	7,0	Width packaging (cm):	11,0
Length (cm) :	13,0	Length packaging (cm) :	16,0
Height (cm) :	4,0	Height packaging (cm) :	5,0
Weight (kg) :	0,130	Weight packaging (kg) :	0,200
SPECIFIC STORAGE CONDITION			
Temperature min	-20	Temperature max	55
Humidity min	5	Humidity max	95



31/08/2022 21:09 - This article is available on our website : www.gsh-med.fr with the following reference SP3110

All characteristics are given for information only and may be changed