

## AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL

Réf. SPG345



# AUTOTENSIO WRIST - WRIST BLOOD PRESSURE MONITOR SPENGLER CORAL



Réf. **SPG345**

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0459	Barcode	3700446041595
Guarantee (month):	24	Customs Code	90189010
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	10,0
Length (cm) :		Length packaging (cm) :	9,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	0,270
SPECIFIC STORAGE CONDITION			
Temperature min	-25	Temperature max	70
Humidity min	5	Humidity max	93



08/12/2022 12:39 - This article is available on our website : [www.gsh-med.fr](http://www.gsh-med.fr) with the following reference SPG345