

DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Réf. SP3110P



DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Réf. SP3110P

FEATURES			
Medical device	No		
Marquage		Barcode	3700446041939
Guarantee (month):	24	Customs Code	90251920
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 03:03 - This article is available on our website : www.gsh-med.fr with the following reference SP3110P

All characteristics are given for information only and may be changed

DISPLAY OF 6 NON-CONTACT THERMOMETER TEMPO PRO SPENGLER



Réf. SP3110P

COMPOSITION

Code	Descriptive	Quantity
SP3110	Non-Contact Thermometer TEMPO PRO Spengler	6