

Réf. DSU231111



Réf. DSU231111

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413	Barcode	3661041008799
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	9,0
Length (cm) :		Length packaging (cm) :	11,0
Height (cm) :		Height packaging (cm) :	14,0
Weight (kg) :		Weight packaging (kg) :	0,250
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:56 - This article is available on our website : www.gsh-med.fr with the following reference DSU231111

All characteristics are given for information only and may be changed

2/2