

Réf. DSU231115



Réf. DSU231115

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413	Barcode	3661041008836
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	25,0
Length (cm) :		Length packaging (cm) :	23,0
Height (cm) :		Height packaging (cm) :	10,0
Weight (kg) :		Weight packaging (kg) :	0,600
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



30/08/2022 21:08 - This article is available on our website : www.gsh-med.fr with the following reference DSU231115

All characteristics are given for information only and may be changed

2/2