

SUNTECH CUFF SINGLE PATIENT SCREW TIGHT 38-50CM

Réf. DSU231059

 SunTech Medical®



SUNTECH CUFF SINGLE PATIENT SCREW TIGHT 38-50CM

Réf. **DSU231059**

FEATURES			
Medical device	Yes	Classe I	
Marquage	CE0413	Barcode	3661041009772
Guarantee (month):	24	Customs Code	
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



31/03/2022 02:58 - This article is available on our website : www.gsh-med.fr with the following reference DSU231059

All characteristics are given for information only and may be changed