

PAIR OF BLACK EARTIPS FOR LAUBRY CLINIC

Réf. 506921



PAIR OF BLACK EARTIPS FOR LAUBRY CLINIC



Réf. 506921

Spare part (x2) for Laubry Clinic® stethoscopes

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage			Barcode	3700446008673
Guarantee (month):	0		Customs Code	90181990
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :			Width packaging (cm):	2,0
Length (cm) :			Length packaging (cm) :	1,5
Height (cm) :			Height packaging (cm) :	1,1
Weight (kg) :			Weight packaging (kg) :	0,005
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



17/11/2022 13:13 - This article is available on our website : www.gsh-med.fr with the following reference 506921

All characteristics are given for information only and may be changed