

Réf. DSU231107

CE
0413



Réf. DSU231107

FEATURES			
Medical device	Yes	Classe I	Accuracy measurement
Marquage	CE0413		Classe I
Guarantee (month):	0	Barcode	3700446014636
Sterile	No	Customs Code	
INFORMATIONS		Main unit	Packaged unit
Width (cm) :		Width packaging (cm):	
Length (cm) :		Length packaging (cm) :	
Height (cm) :		Height packaging (cm) :	
Weight (kg) :		Weight packaging (kg) :	0,378
SPECIFIC STORAGE CONDITION			
Temperature min		Temperature max	(N/A)
Humidity min		Humidity max	(N/A)



31/03/2022 02:02 - This article is available on our website : www.gsh-med.fr with the following reference DSU231107

All characteristics are given for information only and may be changed

2/2