

Réf. DSU250903

CE
0413



Réf. DSU250903

Flexible adult digital SpO2 sensor (90 cm cable)

FEATURES					
Medical device	Yes	Classe I	Accuracy measurement	Classe I	
Marquage	CE0413		Barcode	8698870322289	
Guarantee (month):	0		Customs Code	90181990	
Sterile	No				
INFORMATIONS		Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):		2,0	
Length (cm) :		Length packaging (cm) :		14,0	
Height (cm) :		Height packaging (cm) :		16,0	
Weight (kg) :		Weight packaging (kg) :		0,010	
SPECIFIC STORAGE CONDITION					
Temperature min		(N/A)		Temperature max	(N/A)
Humidity min		(N/A)		Humidity max	(N/A)



31/03/2022 02:02 - This article is available on our website : www.gsh-med.fr with the following reference DSU250903

All characteristics are given for information only and may be changed