

F3000 DISPOSABLE TEMP PROBE COVERS (25 BOXES/TRAY, 20 COVERS/BOX)

 SunTech Medical®

Réf. DSU551467

CE
0413



F3000 DISPOSABLE TEMP PROBE COVERS (25 BOXES/TRAY, 20 COVERS/BOX)

Réf. DSU551467

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	1088452116202
Guarantee (month):	0		Customs Code	
Sterile	No			
INFORMATIONS	Main unit	Packaged unit		
Width (cm) :		Width packaging (cm):	12,5	
Length (cm) :		Length packaging (cm) :	29,0	
Height (cm) :		Height packaging (cm) :	10,0	
Weight (kg) :		Weight packaging (kg) :	0,622	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)	Temperature max	(N/A)	
Humidity min	(N/A)	Humidity max	(N/A)	



31/03/2022 02:03 - This article is available on our website : www.gsh-med.fr with the following reference DSU551467