

Réf. DSU552760



Réf. DSU552760

Spare part for tourniquets

FEATURES				
Medical device	Yes	Classe I	Accuracy measurement	Classe I
Marquage	CE0413		Barcode	3700446033590
Guarantee (month):	3		Customs Code	
Sterile	No			
INFORMATIONS		Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	8,0	
Length (cm) :		Length packaging (cm) :	222,0	
Height (cm) :		Height packaging (cm) :	4,0	
Weight (kg) :		Weight packaging (kg) :	0,270	
SPECIFIC STORAGE CONDITION				
Temperature min	(N/A)		Temperature max	(N/A)
Humidity min	(N/A)		Humidity max	(N/A)



31/03/2022 02:04 - This article is available on our website : www.gsh-med.fr with the following reference DSU552760

All characteristics are given for information only and may be changed