

6X STANDARD BLACK BULBS (LIAN SCENIC AND -S, VL HOPITAL, VL SKIPPER)



Réf. 311130



6X STANDARD BLACK BULBS (LIAN SCENIC AND -S, VL HOPITAL, VL SKIPPER)



Réf. **311130**

Black bulbs

FEATURES			
Medical device	No		
Marquage	Barcode	3700446007492	
Guarantee (month):	60	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	18,0
Length (cm) :		Length packaging (cm) :	12,0
Height (cm) :		Height packaging (cm) :	3,0
Weight (kg) :		Weight packaging (kg) :	0,156
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



23/06/2022 21:04 - This article is available on our website : www.gsh-med.fr with the following reference 311130

All characteristics are given for information only and may be changed