

LIAN CUFF/1 TUBE/NYLON/GREY/ADULT (M)

Réf. 333011



LIAN CUFF/1 TUBE/NYLON/GREY/ADULT (M)



Réf. 333011

Nylon cuff without coupling for Lian® with 60 cm single tubing bladder for Lian® NM, Lian® Métal, Lian® Nano, Lian® Scénic, Ergo Sphyg®, Lian® Color.
Washable at 30-40° C.

FEATURES			
Medical device	Yes	Classe I	
Accuracy measurement	Classe I		
Marquage	Barcode		3700446030377
Guarantee (month):	6	Customs Code	9018901000
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :	Width packaging (cm):		
Length (cm) :	Length packaging (cm) :		
Height (cm) :	Height packaging (cm) :		
Weight (kg) :	Weight packaging (kg) :		0,145
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



25/04/2022 21:06 - This article is available on our website : www.gsh-med.fr with the following reference 333011

All characteristics are given for information only and may be changed