

OSCAR 2 ABP MONITOR - SPHYGMOCOR INSIDE

Réf. DSU230221

SunTech Medical®



OSCAR 2 ABP MONITOR - SPHYGMOCOR INSIDE

Réf. DSU230221

FEATURES			
Medical device	Yes	Classe IIa	
Marquage	CE0413	Barcode	3661040217536
Guarantee (month):	24	Customs Code	90181990
Sterile	No		
INFORMATIONS	Main unit	Packaged unit	
Width (cm) :		Width packaging (cm):	19,0
Length (cm) :		Length packaging (cm) :	31,0
Height (cm) :		Height packaging (cm) :	19,0
Weight (kg) :		Weight packaging (kg) :	1,696
SPECIFIC STORAGE CONDITION			
Temperature min	(N/A)	Temperature max	(N/A)
Humidity min	(N/A)	Humidity max	(N/A)



14/04/2022 21:07 - This article is available on our website : www.gsh-med.fr with the following reference DSU230221

All characteristics are given for information only and may be changed